


Cutting Edge Technology

for the High Definition Operating Room



SDC Ultra

HD Information Management System



Native HD Video Capture

- Record surgical footage in Native HD Resolution
- Perfect for use in teaching, conferences, and presentations
- Unprecedented quality in surgical video



Redesigned Graphic User Interface

- Control the SDC Ultra from the LCD Touch Screen
- Intuitive, easy to use
- Icon Driven user interface



Dual Channel Input Support

- Input/Display two separate video signals
- Capture images and record videos on two separate channels, in synchronized mode, or in PIP format
- On-Screen Picture in Picture to display multiple video signals

SDC Ultra

The Final Piece to the Puzzle



Dual-Rotation Integrated Touchscreen

- Provide flexibility to OR workflow and design
- Dual Rotation Touchscreen for optimal visualization from any angle




Front-Load Patient Scheduling:

- Pre-load patient information for a day's worth of cases
- Easily select patient by convenient drop down menu in the capture screen
- Reduce time and effort between cases for surgical and nursing staff



Easily Customize Annotations, Specialty, & Procedure:

- Annotate images and videos from preloaded list in the album
- Organize annotations and case information by Specialty and Procedure
- Add customized procedures and annotations with one easy step



Full Network Integration:

- Interface with Hospital PACS and FTP Networks
- Save and Retrieve case information from hospital networks/servers
- Auto-Populate patient information from hospital system for accurate medical records



Archive:

- Archive every case with images and videos to the internal 250 GB Hard Drive
- Review, Reprint, and Resave case information from the archive
- Safeguard your hospital - Don't lose another important case